

Texture Analyzer (TA)

Equipment /Instrument Details

Make : Stable Micro-system

Model : TA-HDi, HD-3200

Important features

Texture analyzer (TA-HDi, Stable Micro-system, UK) at ICAR- CARI, Izatnagar operate using Texture Expert Windows Operating System. It is provided with many probes (compressor, shear, puncture etc) and having load cell capacity up to 500 N.



Specifications

Double beam, Load cell 50 Kg, Accuracy 500 PPS

Usages

To determine textural properties of food products

Working Principles

The texture analyzer is work with the principle of compressor, shear, puncture etc and values are determined using Texture Expert software.

Application

For textural analysis of food products

Suggested charges including GST

S.No.	Type of Analysis	Industry	University	ICAR Institutes	Other	Unit
1	Texture	150	100	50	150	Per sample

**The rates are tentative and exclusive of the consumables*

Contact Person
Dr. Ashim Biswas
M. No.: 9463320622
Email ID: biswaslpt@gmail.com

The Demand Draft should be in favour of "The Director, CARI, Bareilly (UP)"
Letter, DD & Samples send to " The Director, CARI, Bareilly (UP)-243122"

Kindly include this note in your requesting letter exactly as it is:

"Content of this report/bill is meant for our information only and we will not use the content of this report for advertisement, evidence, litigation or quote as certificate to third party"